



READYWRAP™
MANUFACTURED BY SOLARIS

Inelastic Medical Binders



ReadyWrap Garments

Thank you for choosing ReadyWrap.

ReadyWrap Professional is a revolutionary, patent-pending step forward in compression garment design. ReadyWrap is constructed with a non-elastic foam-backed fabric which, along with the 50% overlap of straps, mimics traditional bandaging. In addition, once the medical professional has set the garment to the appropriate level of compression, the patient can don and doff the garment while retaining the professionally set compression level.

ReadyWrap Quick-Fit is an economical wrapping system featuring non-elastic straps with a 50% overlap, mimicking traditional, short-stretch bandaging without the hassle.



Table of Contents

ReadyWrap Components & Warranty	4
Washing & Care Instructions	5
Instructions for the Medical Professional	6
ReadyWrap Professional Instructions	10
ReadyWrap Quick-Fit Instructions	14
Checking Compression (Quick-Fit Units)	19
Other Solaris Products	20
Contact Information	24

ReadyWrap Components

- Nylon Laminated Foam
- Latex-free
- Long Lasting Velcro
- Manufactured in the USA

Your Warranty

ReadyWrap Professional garments are warranted to be free of material defects for six months. Exposure to heat over 130° F will damage the garment and void the warranty.

ReadyWrap Quick-Fit garments are warranted to be free of material defects for three months. Solaris will replace any garment found to be defective. Exposure to heat over 130° F will damage the garment and void the warranty.

Washing & Care Instructions

Washing Your ReadyWrap

All ReadyWrap garments should be hand washed using lukewarm water and mild detergent. Using a washing machine or dryer may damage the materials. Lay flat on a towel to dry.

Caring for Your ReadyWrap

To ensure that your ReadyWrap lasts as long as possible, please take note of the following tips:


- Do not don or doff the garment more than advised by your medical professional. The Velcro will wear faster.
- Don't overtighten the garment. The material may stretch.

Professional Setting

Instructions for Setting the Professional Calf Unit (to be Set by the Medical Professional):

Step 1 Open up the rear straps of the unit and make certain the front tabs are all tightly closed (see right).

Step 2 Once the rear straps have been opened, align the center of the front gusset along the patient's tibia, with the Velcro tabs facing out and the smallest straps at the bottom.



Professional Setting

Instructions for Setting the Professional Calf Unit (to be Set by the Medical Professional) continued:

Step 3 Begin fastening the rear straps starting with the bottom strap (the bottom of the garment should be approx. 2" above the ankle).

Step 5 Continue until all rear straps have been set to the desired compression (you may want to have the patient practice donning and doffing once the garment has been set).



Professional Setting

Instructions for Setting the Professional Foot Unit (to be Set by the Medical Professional):

Step 1 Open up the rear straps of the unit and make certain the front tabs are all tightly closed (**see right**).

Step 2 Align the center of the front gusset along the top of the patient's foot, with the straps marked "ankle" (**right**) at the crease of the foot and ankle.



Professional Setting

Instructions for Setting the Professional Foot Unit (to be Set by the Medical Professional) continued:

Step 3 Starting at the toes, begin fastening the rear straps of the unit (see right).

Step 4 Continue until all rear straps have been set to the desired compression (you may want the patient to practice donning and doffing once the garment has been set).



ReadyWrap Professional

Calf Unit Instructions for Donning (after the garment has been professionally fit):

Step 1 Gently slide the ReadyWrap over your foot and onto your calf.

Step 2 Once at the correct height, engage one of the middle straps to hold the unit in place.

Step 3 Starting at the bottom of the garment and working up, engage all the straps.

Step 4 Adjust any straps that did not properly close.



ReadyWrap Professional

Calf Unit Instructions for Doffing (after the garment has been professionally fit):

Step 1 Starting at the top of the garment, release all of the front straps.

Step 2 Gently slide the garment over your ankle and foot to remove.



ReadyWrap Professional

Foot Unit Instructions for Donning (after the garment has been professionally fit):

Step 1 Gently slide the ReadyWrap over your toes and onto your foot.

Step 2 Once properly positioned on the heel, engage the middle two straps to hold the unit in place.

Step 3 Starting at the bottom of the garment and working up, engage all the straps.

Step 4 Adjust any straps that did not properly close.



ReadyWrap Professional

Foot Unit Instructions for Doffing (after the garment has been professionally fit):

Step 1 Starting at the top of the garment release all front straps.

Step 2 Using your hand as a guide, slide the garment over your heel and off the foot.



ReadyWrap Quick-Fit

Calf Unit Instructions for Donning:

Step 1 Align the back spine of the garment with the back of your calf, and the bottom of the garment just above the ankle.

Step 2 Starting with the bottom strap, begin securing each strap to the appropriate level of compression.

Step 3 Continue fastening the remaining straps.



ReadyWrap Quick-Fit

Calf Unit Instructions for Doffing:

Step 1 Starting with the top strap, begin to unfasten each strap in order.

Step 2 As you unfasten each strap, fold the Velcro tabs back onto the strap to avoid tangling.



ReadyWrap Quick-Fit

Foot Unit Instructions for Donning:

Step 1 Align heel in the middle strap heel opening.

Step 2 Fasten the middle strap to the appropriate level of compression.

Step 3 Fasten the strap over the toes to the appropriate level of compression.



ReadyWrap Quick-Fit

Foot Unit Instructions for Donning (continued):

Step 4 Fasten the strap over the ankle to the appropriate level of compression.

Step 5 If needed, adjust straps to ensure the appropriate compression and a comfortable fit.



ReadyWrap Quick-Fit

Foot Units Instructions for Doffing:

Step 1 Remove the ankle strap, folding the Velcro back onto the strap to prevent tangling.

Step 2 Remove the toe strap, folding the Velcro back onto the strap to prevent tangling.

Step 3 Remove the final strap, again folding the Velcro back onto the strap to prevent tangling.



Checking Compression

Checking the Compression of Quick-Fit Units:

Step 1 Donn the garment following the instructions found in this user guide.

Step 2 Beginning at the ankle, slide your index, middle, and ring fingers underneath the strap. If the compression is correct, the strap should only reach the first knuckle of each finger.

Step 3 Continue checking all the straps and adjust any straps which need adjustments until the correct compression is achieved.



Joint Jackets Units

Joint Jackets were developed from the same technology used to treat Lymphedema. Solaris has spent the last four years developing unique products to help athletes recover faster from injuries.

Joint Jackets will:

- Reduce inflammation, swelling, and pain
- Remove lactic acid formation and hematomas
- Allow for customized compression
- Provide support and joint protection
- Increase circulation and massage soft tissues
- Decrease scar adhesions

For more information on Joint Jackets log on to www.solarismed.com



Tribute Therapeutic NightWear Garments

Long-Term Self-Management

Tribute custom-made NightWear garments are essential for a successful lymphedema home management program. Worn nightly, your Tribute garment will offer an effective solution to help regulate fluid fluctuations, decrease tissue turgor, and assist a weakened lymphatic system. Solaris' attention to the details provides you with a garment that eliminates complex bandaging, and is stress free, comfortable to wear, easy to care for, and backed with a one year warranty.



Caresia Bandage Liners

Caresia bandage liners are standard sized bandage liners which simplify complex multi-layered bandaging protocols while enhancing therapeutic outcomes in the active phase of therapy.

Caresia bandage liners are available in full arm, glove, gauntlet, half leg, and foot styles.

To purchase a Caresia, consult your therapist or fitter. Learn more at:
www.solarismed.com



Swell Spots

Swell Spots are professionally designed treatment pads that tuck into garments or under bandaging to address particularly troublesome areas of edema, fill contours, diffuse pressure points, or break up fibrosis.

You can view our entire line of therapeutic products designed specifically to resolve swelling and its complications online at www.solarismed.com



SOLARIS

F R E E D O M T O L I V E

Solaris, Inc.

6737 West Washington Street, Suite 3260

West Allis, WI 53214

Client Hotline: 888-918-9185

T: 888-918-9180

F: 414-918-9189

E: info@solarismed.com

www.solarismed.com